**Mr. Malt®**

A WORLD OF BREWING SOLUTIONS

P.A.B. S.r.l. Via V. Menazzi Moretti, 4 int. 4
33037 - Pasian di Prato (UD) Italy

+39.0432 644279

info@mr-malt.it

www.mr-malt.it

Indice

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www.mr-malt.it

Indice

General	1
Sample preparation	2
Practical application	3
Cleaning	4
Connection to smartphone and tablet.....	5
Instrument check	5

Select language:

EasyDens Digital Density Meter

General

1. What measurement parameters does EasyDens provide?

Starting from version 1.1.1 of the app, EasyDens can provide the following parameters. **Standard measurement parameters**

- Density in g/cm³ and kg/m³, without temperature compensation.
- Specific Gravity (SG), with temperature compensation.

Optional measurement parameters:

- Extract for determining the original and apparent extract in weight (% weight, equivalent to degrees Plato), temperature compensated.
- Alcohol content by volume (% vol), temperature compensated. Used for determining the alcohol content in binary alcohol-water solutions and distillates with low extract content (e.g., liqueurs, gin, vodka, etc.). It is not possible to measure the alcohol content of beer, wine, or similar products due to the high extract content.
- Sugar content (°Brix), with temperature compensation.

In addition to the standard measurement parameters, it is possible to choose the free activation of one of the available optional measurement parameters. The other optional measurement parameters can be activated through an in-app purchase.

2. What are the optional measurement parameters used for?

- Extract (%weight):** measures the extract in beer wort/control of the beer fermentation process with a daily measurement of the apparent extract.
- Alcohol content (%vol):** measures alcohol in binary mixtures of ethanol and water (without added sugars), e.g., distillates, liqueurs, etc.
- Sugar content (°Brix):** measures sugar content in (almost) binary sugar and water mixtures, e.g., fruit juices and various soft drinks. During wine production, a daily measurement of the apparent sugar content allows control of the fermentation process.

3. What comes with EasyDens?

Included in the shipment:

- 1 EasyDens portable hydrometer
- 1 protective case for EasyDens (including 2 o-rings)
- 1 10 mL syringe
- 2 AA batteries
- 2 Luer adapters
- 1 Luer Lock adapter
- 1 drain tube
- 1 illustrated application description
- 1 certificate of conformity
- 1 safety guide

4. Can EasyDens be calibrated with official standards?

EasyDens is not accompanied by any CE certification certificate, so it cannot be calibrated with official standards. In many cases, instruments of this accuracy class are tested for comparison with a laboratory instrument calibrated with official standards or by performing an interlaboratory test.

It is important to note that the inability to calibrate with official standards does not diminish the quality of measurement results. Anton Paar guarantees the specified accuracy and repeatability values.

Sample preparation

1. Do I need to degas the sample before filling the measuring cell?

Yes. To obtain correct and accurate results, the sample must be degassed before measurement.

2. What is the simplest way to degas my sample?

- To degas the sample, choose one of the following methods:
- Shake the sample until no more bubbles form. This can take up to 15 minutes.
- Place the sample in an ultrasonic bath for 5–10 minutes.
- Fill a flask with the sample and seal it with a stopper. Shake the flask and lift the stopper to release the pressure; repeat this procedure until no more pressure is present.
- Fill the supplied syringe halfway with the sample. Then close the syringe tip with your finger and pull back the green plunger. Remove your finger from the tip to release the gas. Push the plunger forward to expel air from the syringe. Repeat these steps several times.

3. Do I need to filter the sample before measurement?

It is not mandatory to filter the sample. However, keep in mind that large particles (>1 mm) can clog the EasyDens inlet and outlet. If the sample contains such particles, it can be filtered using a cotton wool plug. In the case of consecutive measurements of samples from the same batch (e.g., fermentation curve measurements), always treat all samples in the same way.

4. What else should I consider in sample preparation?

Always subject all samples from the same batch to the same preparation treatment. This is the only way to achieve comparable measurement results and avoid potential errors due to different sample handling. Therefore, always use the same procedure for degassing and filtering the sample.

Practical application

1. After switching on EasyDens, I hear a high-pitched tone. Is my instrument defective?

EasyDens works based on the oscillating U-tube method, which determines the density of a liquid from the difference between the oscillations with air in the tube and those with the sample. For this reason, after switching on, the measuring cell is electronically excited to vibrate at its characteristic frequency. The high oscillation frequency can be heard by the human ear as a high-pitched tone. This is not a defect but an indication that the instrument is operating.

3. At what temperature should the sample be for measurement?

The temperature of the sample or the cleaning liquid must be between 0 and 80 °C (32–176 °F). If liquids at higher temperatures are used for filling, the plastic parts may soften and connections may loosen. Also note that automatic temperature compensation only works for sample temperatures between 5 and 30 °C (41–86 °F).

5. For which temperature range does automatic temperature compensation work?

The automatic temperature compensation works for sample temperatures between 5 and 30 °C (41–86 °F).

7. Is it possible to measure with the instrument in a horizontal position?

Yes. Measurements can also be taken with the instrument in a horizontal position. However, to minimize the influence of potential residual gas in the sample, which can cause inaccurate results, it is recommended to take measurements with the instrument in a vertical position. In the vertical position, residual gas in the sample can rise to the upper part of the measuring cell, reducing the risk of inaccurate results.

9. Do I need to drain the sample after measurement?

Yes. After finishing the measurement, the sample must be drained and not returned to the original container. This prevents the risk of contamination.

2. How much sample volume is required?

About 2 ml is needed to fill the measuring cell and perform a correct measurement. In the case of consecutive measurements, it is recommended to fill with 10 ml to completely replace any residual sample from the previous measurement. This amount corresponds to the volume of the supplied syringe.

4. Are the displayed values temperature compensated?

Concentration values calculated from the measured density (e.g., °Brix, % vol alcohol, % weight extract — available from v1.1 onward) are automatically temperature compensated to 20 °C (68 °F). Density values (g/cm³ and kg/m³) are not temperature compensated. For SG (specific gravity), the sample's density compensated to 20 °C is divided by the density of water at 20 °C. A coefficient of 0.0003 g/cm³/K is applied to the sample density's temperature compensation. This coefficient is typical for aqueous solutions.

6. Can I use EasyDens without the protective case?

Yes. EasyDens can also be used without the protective case. However, please note that only with the case installed is the instrument sufficiently protected from splashes, especially during filling. The protective case also makes it easier to use the instrument in an upright position. In the vertical position, residual gas in the sample can rise to the upper part of the measuring cell, reducing the risk of inaccurate results.

8. What should I do after finishing measurements?

After completing measurements, it is recommended to rinse the measuring cell with at least 10 ml of hot distilled or deionized water to avoid contamination or deposits. After rinsing the measuring cell with distilled or deionized water, it is not necessary to dry it since the water will evaporate without leaving any residue. For more detailed cleaning information, refer to section D.

10. Is it possible to perform flow measurements by sending the sample into EasyDens with a peristaltic pump?

Flow measurements are not compliant with the intended use of the instrument. The CE safety report has definitively established, based on the regulations included in section 2.2 of the instruction manual: "Use EasyDens exclusively as a stand-alone unit. Do not combine the instrument with other devices for automated operations."

Cleaning

1. Do I need to clean and dry the instrument after every measurement?

No. It is not necessary to clean and dry the measuring cell between individual measurements. Any residual sample in the measuring cell will be completely replaced by the next sample. EasyDens comes with a 10 ml syringe. It is recommended to use this syringe volume since 10 ml is sufficient to fully replace the previous sample and fill the measuring cell with the new sample. Even if the measuring cell was filled with water, 10 ml of the new sample is enough to avoid dilution. Therefore, there is no risk of incorrect measurement results due to mixing with the previous sample or with water present in the measuring cell.

2. What type of cleaning should I perform after finishing a series of measurements?

After completing a series of measurements, rinse the measuring cell with hot distilled or deionized water. To do this, fill the supplied syringe with hot distilled or deionized water and flush the measuring cell and the drain tube. After rinsing with distilled or deionized water, there is no need to dry it since the water will evaporate without leaving any residue. It is also recommended to clean the measuring cell once a week with a suitable laboratory cleaner (e.g., Mucisol®). Follow the product instructions carefully, especially regarding the detergent concentration. Refer to the EasyDens instruction manual for details.

3. I do not use the instrument regularly. Is there a risk of mold forming in the measuring cell during storage?

If the instrument is not used for long periods, it is advisable to disinfect the measuring cell with high-grade ethanol after rinsing with hot distilled or deionized water.

4. There is a film in the measuring cell that distilled or deionized water does not remove. What should I do?

If visible deposits have formed inside the measuring cell or if calibrations with distilled or deionized water show unacceptable deviations from target values, it is recommended to clean the measuring cell with a suitable laboratory cleaner (e.g., Mucisol®). Follow the product instructions carefully, particularly regarding detergent concentration. Refer to the EasyDens instruction manual.

5. Is it necessary to heat the distilled or deionized water for cleaning?

Especially after measuring wort samples, it is recommended to clean the measuring cell with warm water because it greatly facilitates the removal of sample residues. In this case, it is suggested to first clean the measuring cell with warm tap water, then rinse with hot distilled or deionized water. Ensure no tap water residues remain in the measuring cell, as this can lead to limescale deposits, which may cause incorrect results.

6. Where can I find distilled or deionized water and laboratory cleaner?

Distilled or deionized water is available from any DIY store, pharmacy, or gas station. Suitable laboratory cleaners can be found at laboratory supply stores and from suppliers of brewing equipment.

7. Is distilled or deionized water suitable for cleaning and calibrating the instrument?

Generally, the quality of water from a DIY store or a gas station is sufficient, provided the water is explicitly labeled as distilled and deionized.

Connection to smartphone and tablet

1. Which types of smartphones and tablets are compatible with EasyDens?

EasyDens is compatible with any smartphone or tablet equipped with Bluetooth LE (also known as Bluetooth Low Energy, BLE, Bluetooth Smart, or Bluetooth Smart Ready). In addition, the device's operating system must be Android 4.3 (or higher) or iOS version 7.0 (or higher).

3. How do I connect my mobile device with EasyDens?

1. Download the free Anton Paar EasyDens app on your smartphone or tablet from Google Play Store or Apple App Store.
2. Enable Bluetooth on your smartphone or tablet.
3. Switch on EasyDens by pressing the on/off button once.
4. Launch the EasyDens app on your smartphone or tablet.
5. The smartphone or tablet will automatically connect to EasyDens.

Instrument check

1. How can I verify that EasyDens is working?

EasyDens is compatible with any smartphone or tablet equipped with Bluetooth LE (also known as Bluetooth Low Energy, BLE, Bluetooth Smart, or Bluetooth Smart Ready). In addition, the device's operating system must be Android 4.3 (or higher) or iOS version 7.0 (or higher).

3. How can I check if EasyDens measurement values are correct?

It is sufficient to compare the density values provided by the instrument for distilled or deionized water with the correct values (at the respective temperature). The correct values can be found in Appendix B of the EasyDens instruction manual. There is also a detailed explanation of the control measurements and the required steps.

2. Is an active connection required to use the app?

The EasyDens app does not require an active Internet connection for normal use to control EasyDens. However, an active Internet connection is required to download the app from Google Play Store or Apple App Store.

4. My smartphone/tablet does not recognize EasyDens

Follow the steps below to resolve connectivity issues between EasyDens and your smartphone or tablet:

1. Check your smartphone or tablet specifications to ensure they meet the basic requirements for communication with EasyDens (see question E.1).
 2. If the connection between the device and EasyDens switched on is not automatically established when launching the app with Bluetooth enabled, return to the app's start screen and tap the "Measure" button again. The smartphone or tablet will restart the search for EasyDens.
 3. If a connection is still not established, check the Bluetooth visibility settings on your smartphone or tablet. Make sure the device is visible to other devices and allows connections.
 4. For proper connection to Bluetooth LE devices, your smartphone's approximate location is required. Check the location settings on your smartphone and make sure the smartphone is not completely switched off (the procedure also works with energy-saving mode active).
 5. If the necessary app permissions (e.g., image access) were refused during the installation procedure, the app will not work correctly. In this case, uninstall the app and reinstall it, granting the required permissions.
- If all the steps above do not allow establishing the connection between the device and EasyDens, send a log file to Anton Paar by email using the relevant function.

2. Do I need to check EasyDens regularly?

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4. There is a small black spot on my measuring cell?

The small black spot on the curve of the measuring cell is the NTC temperature sensor. Since the density of a liquid depends on temperature, the EasyDens temperature sensor is placed directly inside the measuring cell to guarantee accurate temperature compensation.

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